



Product information

Cryo Box

Catalog #: SP91192, SP91292
Size: 10 pc (SP91192), 40 pc (SP91292)

1. Description

Cryo Boxes are useful for storing vials in ultra low temperature environments without distortion of storage boxes.



2. Features

- For storing of 1.5 ml microtubes or cryogenic vials
- Designed for enhanced drainage of liquid nitrogen (SP91192).
- Alphanumeric indications for sample location.
- Temperature range: -196 ~ +121°C (SP91192).
- 3 colors: Blue, Pink, Green (SP91192).
- Cardboard Box: Moisture repellent coating (SP91292). Designed to fit both 1.2 ml / 1.8 ml Cryovials.

Cat. #	Material	Feature	Description	Packaging
SP91192	PC	3 colors	9 x 9 (81 Holes)	10/case (1/sleeve)
SP91292	Cardboard	-	9 x 9 (81 holes)	40/case

3. Dimensions

Cat. #	Length (mm)	Width (mm)	Height (mm)	Material
SP91192	133.5	133.5	35.2	Polycarbonate (PC)
SP91292 CAP	132.5	132.5	35	Cardboard
SP91292 BODY	127.4	127.4	50	Cardboard



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
NOT INTENDED FOR HUMAN OR ANIMAL USE.