



CERTIFICATE OF ANALYSIS

Product	0.85% Sodium Chloride-Phosphate Buffer pH 7.3 (770L)
Product Code	A212107-0770
Lot No	
MFG Date	
Manufacturing Site	Canada
Shelf Life	18 Months

Test Items	Specifications	Results
Bioburden	≤ 3 CFU/30 mL	
Endotoxin	≤ 1.00 EU/mL	
pH	7.2 - 7.4	
Conductivity	18.1 – 19.6 mS/cm	

Composition: Per (770 L)

770kg Water for Injection

8100.4g Sodium Phosphate Dibasic, Heptahydrate (CAS # 7782-85-6)

6545g Sodium Chloride (CAS # 7647-14-5)

1062.6g Sodium Phosphate Monobasic, Monohydrate (CAS # 10049-21-5)

Storage: Room Temperature (15 – 30 °C).

Filtration Method: Filtered through a 0.22µm filter.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: April 26, 2024