



CERTIFICATE OF ANALYSIS

Product	0.05M PBS, 0.85% Sodium Chloride, pH 7.3 (500L)
Product Code	A212109-0001
Lot No.	
MFG Date	
Manufacturing Site	Canada
Shelf Life	18 Months


Test Items	Specifications	Results
Endotoxin	≤ 0.25 EU/mL	
pH	7.2 - 7.4	
Osmolality	361 – 399 mOsm	

Composition: Per (500 L)

- 500kg Water for Injection
- 5260g Sodium Phosphate Dibasic, Heptahydrate (CAS # 7782-85-6)
- 690.0g Sodium Phosphate Monobasic, Monohydrate (CAS # 10049-21-5)
- 4250g Sodium Chloride (CAS# 7647-14-5)

Storage: Room Temperature (15 – 30 °C).
Filtration Method: Filtered through a 0.22µm filter.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: April 26, 2024