



CERTIFICATE OF ANALYSIS

Product	50mM Tris Buffer (150 KG)
Product Code	A212110-0001
Lot No	
MFG Date	
Manufacturing Site	Canada
Shelf Life	18 Months

Test Items	Specifications	Results
Bioburden	≤ 3 CFU/30 mL	
Endotoxin	≤ 0.25 EU/mL	
pH	Report	
Osmolality	Report	

Composition: Per (150 KG)

- 141 kg Water for Injection
- 908.5 g TRIS (CAS # 77-86-1)
- 620 g Hydrochloric Acid (CAS # 7647-01-0)

Storage: Room Temperature (15 – 30 °C).
Filtration Method: Filtered through a 0.22µm filter.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: April 26, 2024