

## Product information

QF 24 V4  
V1 Sept 2024

# PSB Powder

**Catalog #:** PD0100  
**Size:** 1 PK  
**Storage:** 18~25°

### Product Description:

This PBS formulation is recommended for use in immunohistology protocols.

Each package makes 1 L of 10x concentrated PBS or 10 liters of 1x solution. Each 1x PBS solution contains 137 mM Sodium Chloride, 2.7mM Potassium Chloride and 10mM Phosphate Buffer, pH 7.3-7.5(1.0%, Water) at 25°C.

### Storage:

Store at 18~25°C.

### Precautions and Disclaimer:

For Laboratory Use Only. Not for drug, household or other uses.

**MUST USE ENTIRE CONTENTS FROM THE CONTAINER PER PREPARATION. The powders are mixed however not blended.**



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.  
NOT INTENDED FOR HUMAN OR ANIMAL USE.