



## CERTIFICATE OF ANALYSIS

---

Product	10x PBST Solution (10X Phosphate Buffered Saline, 0.5% Tween-20 <sup>TM</sup> of ICI America), pH 7.4
Grade	Cell culture
Product Code	UPD8119
Lot No	
Re-Assay Date	

---


---

Test Items	Specifications	Results
Appearance	Liquid	
pH	7.3~7.5	
Sterility	Pass	
DNase	None Detected	
RNase	None Detected	
Protease	None Detected	

---

Storage: 18-25°C.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature:   
Title: General Manager  
Date: November 11, 2024