



Product Information

BT999-SN

(5.0ml Tip, Natural, Sterile, 50/Pack)


Features	Specifications
Item Number	BT999-SN. SIZE .1PK
Color	Natural
Material	PP
Tip (Volume capacity)	5.0 ml
Sterile or Non-Sterile	Sterile
Tips per Pack	50
Dimensions	Please see picture on Page 2
Storage	Room Temperature



Picture: 1



MAT'L: PP	GEL.TOL.: ±0.05	PROG:	DRW.NO.: BT620	APPR. BY:
DATE: 2012/12/29	UNIT: MM	SCALE: 1/1	DRAWN BY: lu_yh	



Picture: 2