

# TIPS SIZE SPEC

Size													
Feature	Non-Graduated	Graduated	Extra Long Graduated	Non-Graduated	Graduated		Non-Graduated	Graduated		Graduated	Extra Long Graduated	Graduated	Non-Graduated
Model	A	B	C	D	E		F	G		H	I	J	K
Standard Tips Cat. No.	F600215 F601215 F602215	F600216 F601216 F602216	F602218	F600225 F601225 F602225	\	F600224 F601224 F603224 F604224	F600227 F601227 F603227 F604227	\	F600330 F601330 F602330	F600222 F601222 F602222 F603222	F604231	F600500 F601500	F501010
Filter Tips Cat. No.	F610215 F611215 F612215	F610216 F611216 F612216	\	F610225 F611225 F612225	F610224 F611224 F612224	\	\	F610226 F611226 F612226	\	F610222 F611222 F612222	F612231	\	\
Volume	10 $\mu$ l			20 $\mu$ l	100 $\mu$ l	200 $\mu$ l		300 $\mu$ l	1 ml		5 ml	10 ml	

Noted: Due to the addition of filter mambrance, some filter tips have a reduced capacity compared to standard tips of the same model. For example, 200  $\mu$ l and 300  $\mu$ l standard tips will become 100  $\mu$ l and 200  $\mu$ l filter tips after adding filter.

